


COVID-19 | ANTIGEN RAPID TEST

Test for self-testing by laypersons

Art.No. 01760



Scope of delivery.

1 pack contains:

5 individually packaged rapid tests.

- 5 Test Cassettes.
- 5 Tubes.
- 5 Swabs.
- 1 Instruction for use.
- 5 Extraction Reagents.

5
tests

Test for self-testing by laypersons

Manufacturer: Hangzhou Clongene Biotech Co., Ltd.





File number of the special approval of the BfArM	Manufacturer	Applicant	Testname
5640-S-168/21	Hangzhou Clongene Biotech Co., Ltd.	Hangzhou Clongene Biotech Co., Ltd.	COVID-19 Antigen Rapid Test

https://www.bfarm.de/DE/Medizinprodukte/Antigentests/_node.html

Product details.

A lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasal swab from individuals aged 7 years or older.

The test is particularly characterized by its simple and safe application. To perform the test, a swab is taken from both nostrils and evaluated using the test cassette. After only 15-20 minutes, the test shows the result.

-  Simple sample extraction.
-  No additional equipment required.
-  Results easy to read.
-  Fast results.

Subject to alterations and errors.