



Quality  
Made in  
Germany.

# Respiratory Mask „CareOne“ FFP2 NR<sup>01769</sup>

**CERTIFIED**

In accordance with Regulation (EU) 2016/425  
and the European EN 149:2001 + A1:2009 standard.


0757 **CE**

## Our quality promise.

- **Machine, material, production, testing institute - 100% Made in Germany.**
- High-quality processing in accordance with European standards.
- Highest standards with regard to environmentally-friendly and hygienic production processes, packaging and storage.
- Uncompromising and consistent quality checking for the most exacting requirements and security

## Your benefits.

- The use of high-quality German filter fabrics means that there is no need for redundant additional filter layers. This makes breathing and the long-term wearing of the mask much more pleasant than with the majority of other models.
- Safe and effective protection for the wearer against particles, liquids and solid aerosols.
- High-quality mask with 4 protective layers. Free from glass fibres, latex and odours.
- Extra skin-friendly non-woven fabric with very low breathing resistance, a higher level of comfort and an exceptional fit.
- **Filtration efficiency (according to standard)  $\geq 94\%$**



Extra skin-friendly non-woven material

Spunbonded non-woven material

Meltblown non-woven material

Spunbonded non-woven material



4-layer  
protection



Aerosol  
protection



Breathable &  
odour-neutral



Elasticated  
ear loops



Flexible  
nosepiece



[elastocare.com](https://elastocare.com)